

DARK CIRCLES AND "TEARS VALLEY" THE ADDED VALUE OF HX

Dark circles are very common in women and men, they depend on the age and the intrinsic characteristics of the person in fact, even in youth dark circles are accompanied by a dark color on the skin.

The so-called "Circles" may over time certainly accentuate or be added to "Bags" and "tears valley" (between cheekbones and nose).

The most immediate solution is definitely blepharoplasty surgery, but this choice is not for everyone considering the costs and post-operative recovery time.

Recently, the less expensive alternative is certainly represented by injectable cross-linked (or not) biomaterials based on hyaluronic acid (for its safety and easy manageability) and calcium hydroxyapatite preferred (preferred for density, white color, and for the "biological correction" induced by the hydroxyapatite micro particles).


Peptidyal HX® is the ideal alternative solution for treating such imperfections. It includes in its formulation, all the positive aspects of hyaluronic acid and hydroxyapatite, avoiding their crosslinking that can sometimes "stimulate" undesired reactions even if protein residues equal to zero.

SUPERIOR PURITY OF RAW MATERIALS

It is known that the implant under eyelid is very sensitive, it is important therefore not to exceed the amount of injected product in order to avoid the appearance of edema, which despite being reversible and definitely temporary, they can persist for several days and be very annoying to the patient.

To avoid this problem, **Peptidyal HX®** has been designed with an extremely "natural" approach using a particular **hyaluronic acid obtained from a non-pathogenic bacterial strain completely devoid of protein residues.**

This produces a product whose safety is superior to the standard.



Recommended Treatments:

- After few treatments **Peptidyal HX®** formulation provides compactness and tonicity promoting hydration and densification of face and body skin.
- 3 biweekly sessions are recommended to be alternated with other products of the Peptidyal range (**Peptidyal®** and **Peptidyal 2%®**) depending on the conditions.
- It is suggested, depending on the areas to be treated, to be used with a series of microbump infiltrations or reticular technique (multiple injections in horizontal and vertical direction to form a grid).
- **Before and after laser treatment**

Each pack contains 2 pre-filled syringes of 2.5 ml each

Composition: Sodium Hyaluronate 18 mg / ml Hydroxyapatite CaHA 0.01% Glycine and L-Proline



Exclusive World Distributor

Aerazen Lab s.r.l.
Viale Majno 34 - 20129 Milano - Italy
info: +39 392 97 46 341 - info@aerazenlab.it
www.aerazenlab.it

Exclusively Manufactured for



PEPTIDYAL | HX

EXCLUSIVE BIO REGENERATIVE FORMULA

A NEW CONCEPT OF BIOSTIMULATION AND BIOREVITALIZATION



Peptidyal HX® can boast the currently available highest concentration of Hyaluronic Acid (18mg/ml), combined with aminoacids.

The ingredients that make up this innovative **injectable medical device (skin booster)**, consisting of a transparent viscoelastic hydrogel, perform an innovative synergy ideal for both face and body treatments.

Non-animal Hyaluronic Acid, is obtained through an innovative manufacturing process that ensures a unique and **extraordinary purity due to the non-use of chemical solvents** as normally done for extracting synthesis hyaluronic acid.

The Hyaluronic Acid of **Peptidyal HX®** has a specific medium molecular weight obtained after computing complex stoichiometric coefficients which allows us to obtain:

- a superficial action, by binding to the stratum corneum and ensuring exceptional hydration and filmogenic action with a tensor and protective effect
- a deep action by stimulating proteins production responsible for compactness and elasticity in the deepest layers of the skin thus resulting in wrinkles and skin depression reduction.

Proline and glycine amino acids in combination with hyaluronic acid help the production of new collagen considering that, as the age progresses, fibroblasts require an increasing amount of amino acids.

Peptidyal HX® contains also **0,01% of calcium hydroxyapatite**: the presence of this important innovative ingredient in a biostimulative product is justified by the fact that **hydroxyapatite microspheres allow fibroblasts to produce new collagen and the Hyaluronic acid degradation is more gradual** due to the presence of CaHa.

CLEAR AND VISIBLE ANTI AGING RESULTS SINCE THE FIRST TREATMENTS

Gender: Female
Age: 36 years
Affected area: Lip contour
Protocol: 2 every 3 weeks
Product (s): Peptidyal HX® + Peptidyal®
Infiltrations: microbump



Gender: Female
Age: 67 years
Affected area: Wrinkled Eyes + Cheeks
Protocol: 3 on a bi-weekly basis
Product (s): Peptidyal HX®
Infiltrations: microbump



Gender: Female
Age: 47 years
Affected area: Cheeks
Protocol: 3 on a bi-weekly basis
Product (s): Peptidyal HX® + Peptidyal 2%®
Infiltrations: microbump



THE 4 FEATURES THAT MAKE PEPTIDYAL HX® A UNIQUE PRODUCT

1

exceptionally pure hyaluronic acid

2

very high concentration (18mg/ml) of HA

3

calcium hydroxyapatite (CaHa) in microspheres

4

proline + glycine + hyaluronic acid = new collagen